

NASA TECH BRIEF

Ames Research Center




NASA Tech Briefs announce new technology derived from the U.S. space program. They are issued to encourage commercial application. Tech Briefs are available on a subscription basis from the National Technical Information Service, Springfield, Virginia 22151. Requests for individual copies or questions relating to the Tech Brief program may be directed to the Technology Utilization Office, NASA, Code KT, Washington, D.C. 20546.

Mass Separator for Low Velocity Ions

The problem:

To provide a broad beam of protons free of other charged species in order to simulate the solar wind in experiments that are designed to determine its effect on materials.



The solution:

Purify the effluent beam from an RF discharge or other source operating with hydrogen by separating the protons from the H_2^+ and H_3^+ ions with the aid of a magnetic lens.

How it's done:

The separator is based on the fact that the focal length of a magnetic lens depends on the momentum of the charged particle that traverses it.

Referring to the figure, it is evident that only particles of a particular mass will be focused through the aperture. Other masses will be focused at different axial positions and thus will form a larger image at the plane of the aperture. When properly designed, the magnetic lens system will permit essentially all

ions of the desired mass to pass through the aperture. Only a small fraction of the others will be transmitted along the axis of the beam. These unwanted ions plus charge-exchange neutrals and photons from the source can be more or less completely eliminated by a small opaque stop placed on axis. Further purification of the beam could be achieved by adding another lens and aperture in series with the first.

The large $\Delta m/m$ associated with low-mass ions aids in making this simple system highly efficient for the purification of proton beams. The magnetic lens system is small and lightweight and, moreover, it requires only a single DC power supply; it can operate at high voltage and thus be used to separate the ions before acceleration. Transmission approaches 100%, compared to less than 10% for the RF separator commonly used for this application.

Notes:

1. This type of mass separator should be generally useful for separating all types of low-energy (<10 keV) light ions.
2. The following documentation may be obtained from:

National Technical Information Service
Springfield, Virginia 22151
Single document price \$3.00
(or microfiche \$0.95)

Reference: NASA CR-73443 (N70-40901),
Solar Wind Simulation Techniques

3. No additional documentation is available. Specific questions, however, may be directed to:
Technology Utilization Officer
Ames Research Center
Moffett Field, California 94035
Reference: B72-10123

(continued overleaf)

Patent status:

No patent action is contemplated by NASA.

Source: Harry J. King of
Hughes Aircraft Company
under contract to
Ames Research Center
(ARC-10375)